

### Battery Operated Cannulated Drill System

It is used for various types of drill operations. Also suitable for K-wire Surgery & Drill/Reaming Operations

Consists of the following:

- Drill Handpiece - 1 Pc.
- Battery - 2 Pcs.
- Charger with cord - 1 Pc.
- Jacob's Chuck with Key - 1 Pc.
- Sterilizing Rings - 2 Pcs.
- Carry Case - 1 Pc.



Battery Operated

### SPECIFICATIONS

#### ★ Handpiece

Speed	1000 rpm
Noise	65db
Sterilize	135 °C
Torque	3.5 NM

Autoclavable

Packed in Aluminium Carrying Case

#### ★ Battery

Type	Ni-MH
Voltage	14.4 V
Charging time	2 Hrs
Discharge time	45-60 Min

#### ★ Charger

Input voltage	100V~240V
Input power	50Hz~60Hz
Input current	0.5A (max)
Output voltage	7.2V~14.4V
Output current	1.0A (max)

### ACCESSORIES



AO Reamer Adapter



Quick Coupling Drill Bit Attachment



Wire/Pin Driver Attachment



Battery



Charger

Better Battery life

Faster Speed

High Power

Better Grip


Ergonomical Shape

### Battery Operated Oscillating Saw System

It is used for various types of Joints Replacement Surgeries.

*Consists of the following:*

- Saw Handpiece - 1 Pc.
- Battery - 2 Pcs.
- Charger with cord - 1 Pc.
- Blades - 2 Pcs.
- Sterilizing Rings - 2 Pcs.
- Carry Case - 1 Pc.



#### SPECIFICATIONS

##### ★ Handpiece

Speed	0-20000 osc/m
Noise	70db
Sterilize	135 °C

Autoclavable

##### ★ Battery

Type	Ni-MH
Voltage	14.4 V
Charging time	2 Hrs
Discharge time	45-60 Min

Packed in Aluminium Carrying Case

##### ★ Charger

Input voltage	100V~240V
Input power	50Hz~60Hz
Input current	0.5A (max)
Output voltage	7.2V~14.4V
Output current	1.0A (max)

#### ACCESSORIES



Notch Blade



broad Blade

Better Battery life

Faster Speed

High Power

Better Grip

Ergonomical Shape