

RZ Medizintechnik GmbH
Unter Hasslen 20 | D-78532 Tuttlingen
Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50
sales@rz-medizintechnik.com
www.rz-medizintechnik.com



RZ
Medizintechnik
GmbH

Progress in Surgery

Laparoscopic LIVER RETRACTOR



news
06.2015

NATHANSON RETRACTOR SYSTEM

The Nathanson system allows the surgeon to hold the entire liver. It is used primarily in Gastroesophageal Reflux Disease procedures, but can also be used in a wide range of laparoscopic surgeries. With the additional accessories, it extends its supporting features for open surgeries, even for larger or bariatric patients. The retractor system is available as a laparoscopy set, each open surgery item can be ordered separately.

NATHANSON LIVER RETRACTOR SET

- Laparoscopic liver holding system
- Quick lock
- Expandable
- Very adjustable
- Bariatric patients hooks



170-801-500 complete set

Order Number	Description
170-801-501	Holding Device, 90°, 360° rotatable
170-801-503	Table Clamp, heavy duty
170-801-506	Standard Instrument Clamp
170-801-507	Scope Holder, insulated
170-801-510	Nathanson Hook, small
170-801-511	Nathanson Hook, medium
170-801-512	Nathanson Hook, large
170-801-515	Nathanson Tray, with lid


Order Number	Description
170-801-503	Table Clamp, heavy duty



Order Number	Description
170-801-507	Scope holder, insulated



Order Number	Description
170-801-506	Standard Instrument Clamp





Quick lock mechanism

Hook, Ø 6.5mm

Order Number	Description
170-801-510	Small
170-801-511	Medium
170-801-512	Large



Order Number	Description
170-801-501	Holding Device, 90°, 360° rotatable



Order Number	Description
170-801-502	Double Holding Device, 90°, 360° rotatable*

*Item is not included in the set 170-801-500

ACCESSORIES FOR OPEN SURGERY




Order Number	Description
170-801-504	Universal Retractor Holder



Order Number	Description
170-801-505	Frame Clamp for Walter frames



Order Number	Description
170-801-508	Ratchet mechanism



Complete Nathanson
Liver Retractor Set is shown
170-801-500

Order Number	Description
170-801-515	Tray, with lid, 540 x 253 x 197mm

© 06/2015 Copyright by RZ Medizintechnik GmbH · Unter Hasslen 20 · D-78532 Tuttlingen. All rights reserved.
Modifications of technical specification or design are subject to change without prior notice.

TriMed Technologies

Customise Surgical Solution & Instruments

4/E-1, Second Floor, Nuchana Tower,
Alaknanda Shopping Complex, Alaknanda, New Delhi - 110019 (INDIA)
Ph:- +91-11-41051188, Mobile:-+91-8459130431
Web: www.trimed.co.in
Email: trimedtechnologies@gmail.com,info@trimed.co.in

