



IriMed

**Ophthalmology
Instruments**



**ISO & CE Certification
For More Reliability...**



Quality Policy

"We are committed to provide a better quality product which suits the customer requirements & keep on improving resource utilization to best of our efforts"

INDEX**Titanium****INSTRUMENTS****PAGES**

Forceps	1,2,3,4,5
Needle Holder	5,6
Scissors	6
Speculum	6,7
Folder	7
Inserter	8
Injector	8
Prechopper	8
Repositor	8
Spatula	8,9,10
Manipulators	8,9
Chopper	9
Divider	9,10
Irrigation Chopper/Y Rotator	10
Bimanual	10
Retractor	11
Vectis	11
Breaker	11
Micro Surgical Sets	12

Stainless Steel

Rotators / Phaco Choppers / Pre Choppers	13,14
Folder / Inserter	15
Manual Phaco / Nucleus Dividers	16
Sapphire Ophthalmic Blades	17
Phaco Knives & Keratomes	17
Speculums	18,19,20
Scleral Depressors / Retractors	20,21
Retractors / Caliper Fixation Rings	21
Trephines / R.K. Instruments	22
Curettes Polisher / IOL Hooks	23
Manipulators	24
Spuds / Loops / Cystotomes	25
Spoon / Spatula / Repositors & Knives	26
Muscles Hooks / Cautery / Expressors	27

Titanium

Muscles Hooks / Cautery / Expressors 27
 Blades Breakers / Dissectors & Bone Punches / Nibblers.... 28
 Bone Punch Nibblers / Lacrimal Instruments 29
 Probes / Dilators & SQ Handle Forceps / Scissors 30
 Clamps / Forceps 31,32,33,34,35,36,37,38
 Scissors 38,39,40,41
 Needle Holder 41,42,43
 Cannulas 43, 44, 45, 46, 47
 Handles / Trays 47
 Lasik / Lasek 48,49,50

Stainless Steel


SET-1 Intra/extra Capsular Cataract with IOL Implataion Set 51
 SET-2 Non-Phaco Set (Small Incision) 52
 SET-3 Non-Phaco Set (Small Incision) with Sapplire surgical blade.52
 SET-4 Phaco Emulsification Set 53
 SET-5 Lacrymal / D.C.R. Set 54
 SET-6 Muscle Surgery Set 55
 SET-7 Glaucoma Set 56
 SET-8 Corneal Transplant Set 57
 SET-9 Enucleation Set 58
 SET-10 Vitrectomy Set 59
 SET-11 & SET - 12 Lasik & Lasek Set 60

New


New Instruments 60,61,62
 Equipments 64,65,66,67,6,69,70,71,72

Forceps


Forceps tying straight with platform, overall length 87mm
MF-30



Forceps tying straight Castroviejo with platform overall length 109mm
MF-30A



Forceps tying straight kirby model with platform overall length 103mm
MF-33



Forceps tying curved with platform, overall length 87mm
MF-36



Forceps tying curved with platform, Castroviejo overall length 109mm
MF-36A



Forceps tying curved kirby model with platform overall length 103mm
MF-39



Forceps tying colibri with platform overall length 83mm
MF-42



Forceps McPherson angled 10mm, with platform overall length 87mm
MF-48



Forceps McPherson angled 12mm Castroviejo with platform, overall length 109mm
MF-48A



Forceps McPherson angled kirby model with platform overall length 103mm
MF-51



Forceps McPherson angled 12mm with platform overall length 85mm
MF-54





Forceps

Forceps
capsulorhexis straight
overall length 85mm
MF-60



Utrata Forceps
capsulorhexis Castroviejo
overall length 109mm
MF-60A



Utrata Forceps capsulorhexis
kirby round handle
overall length 103mm
MF-63



Utrata Forceps
capsulorhexis curved shaft
kirby round handle
overall length 103mm
MF-63A



Utrata Forceps
capsulorhexis curved shaft
overall length 85mm
MF-66



Utrata Forceps
capsulorhexis curved shaft
castroviejo handle
overall length 109mm
MF-69



Capsulorhexis Forcep
23 gauge Curved Shaft
with micro grasping tip
MF-69A



Corneal Forceps toothed 0.12mm/0.2mm,
1x2 teeth with tying platform straight
overall length 87mm
MF-72



Corneal Forceps toothed
0.12mm/0.2mm, 1x2 teeth
with tying platform straight
castroviejo handle
overall length 107mm
MF-72A



Corneal Forceps toothed
0.12mm/0.2mm, 1x2 teeth
with tying platform kirby
round handle, overall length 103mm
MF-75



Corneal Forceps toothed 0.12/0.2mm,
1x2 teeth with tying platform curved
overall length 87mm
MF-78



Corneal Forceps toothed
0.12mm, 1x2 teeth with tying
platform/curved kirby round
handle, overall length 107mm
MF-78A



Forceps

Corneal Forceps colibri
0.12mm, 1x2 teeth with tying
platform kirby round handle
overall length 103mm
MF-81



Colibri Forceps 0.12mm,
1x2 teeth with tying platform
overall length 83mm
MF-84



Polock Corneal Suturing forceps
very delicate double
1x2 teeth 0.12mm
overall length 83mm
MF-84A



Lims Forceps 0.12mm/0.2mm,
1x2 teeth with tying platform
overall length 93mm
MF-87



Lims Forceps pierse
type with tying platform
overall length 87mm
MF-87A



Lims Forceps 0.12mm/0.2mm,
1x2 teeth with tying platform
overall length 87mm
MF-87B



Forceps notched straight
pierce with tying platform
overall length 87mm
MF-90



Forceps notched curved
pierce with tying platform
overall length 87mm
MF-96



Forceps notched colibri
pierce with tying platform
overall length 83mm
MF-102



Forceps pierse type with
tying platform kirby round
handle, overall length 103mm
MF-105



Forceps capsule tag
1.5mm platform at tips
overall length 85mm
MF-108





Forceps

Forcep eye dressing
serrated delicate,
overall length 100mm
MF-111



Forcep fixation 1x2 teeth
overall length 90mm
MF-114



Forcep
superior rectus toothed
overall length 115mm
MF-117



Inomura
Capsulorhexis Forceps
overall length 115mm
MF-120



Chalozion
Hunt Forceps
overall length 95mm
MF-123



Chalozion Forceps
Oval Shape
overall length 95mm
MF-126



Jewellers Forceps
overall length 100mm
MF-129 Straight Fine
MF-132 Straight with platform



Faulkner Soft IOL
Holding Forceps
overall length 109mm
MF-135



Lens holding
forceps 'T' type
overall length 110mm
ML-03



Lens holding forceps
McPherson type
overall length 103mm
ML-09



Lens holding
forceps 'L' type
overall length 110mm
ML-15



Forceps / Needle Holder


Lens holding
forceps kratz
smooth curved jaws
overall length 110mm
ML-21




Lens holding
forceps Daljit model
overall length 85mm
ML-24



Needle holder barraquer
straight without lock
overall length 113mm
MN-300



Needle holder barraquer
straight with lock
overall length 113mm
MN-306



Needle holder barraquer
curved without lock,
overall length 113mm
MN-312



Needle Holder
Baraquer W/o Lock
Accurate Curved
overall length 113mm
MN-315



Needle holder barraquer
curved with lock,
overall length 113mm
MN-318



Needle holder
barraquer long model
overall length 135mm
MN-321



Needle holder
castroviejo curved
delicate without lock
overall length 113mm
MN-336





Needle Holder / Scissors / Speculum

Needle holder castroviejo
 long curved
 overall length 120mm
MN-342 with lock
MN-345 without lock



Needle holder castroviejo
 long straight
 overall length 120mm
MN-348 with lock
MN-351 without lock



Scissor corneal
 universal curved
 overall length 100mm
MS-240



Scissor corneal straight
 overall length 100mm
MS-243



Scissor vannas curved
 short blade/ long blade
 overall length 95mm
MS-246



Scissor vannas straight
 overall length 95mm
MS-249




Scissor vannas angled
 upward capsulotomy
 overall length 95mm
MS-252



Speculum, adjustable
 overall length 82mm
MS-255 with 15mm solid blade for adults.
MS-255A with 11mm solid blade for infants & child



Speculam adjustable spread
 for phaco / lasik
 overall size 82mm
MS-258



Speculum / Folder / Inserter


Lieberman Speculum adjustable
open wire blade
MS-261



Speculum adjustable
open wire blade S.S.
MS-264



Kratz-Barraquer Wire Speculum
with open blades for unobstructed
access to the cataract incision adult / child
MS-267



Phaco silicon
soft lens folder
overall length 108mm
MP-21



Phaco silicon soft lens
insertor with lock
overall length 115mm
MP-27



Phaco acrylic lens folder
overall length 110mm
MP-30



Phaco acrylic lens insertor
overall length 108mm
MP-33



Phaco acrylic soft
lens folder box type
overall length 116mm
MP-36



Phaco acrylic soft lens
insertor to be used with
MP36 lens folder
overall length 112mm
MP-39



Folder Phaco acrylic soft IOL's
folding forceps with smooth platform
overall length 111mm
MP-42



Insertor for folding & holding the
shaft IOL's for insertion
overall length 100mm
MP-45





Injector / Prechopper / Repositor / Spetula / Manipulator

Injector for acrylic soft lenses
overall length 121mm
MP-48



Pre-Chopper, Angled Round Handle
overall length 121mm
MP-54



Pre-Chopper, Curved Round Handle
overall length 121mm
MP-57



Akakoshi Pre-Chopper for cracking the soft or hard nucleus
overall length 121mm
MP-60



Repositor single ended malleable
overall length 123mm
MR-204



Spatula single ended
overall length 123mm
MR-207



Medelec repositor curved/angled, double ended,
overall length 145mm
MR-207A



Sinsky hook extra fine single ended
overall length 123mm
MR-210



Sinsky hook cum spatula double ended
overall length 145mm
MR-213



No hole IOL manipulator
overall length 115mm
MR-216



IOL manipulator ball ended
overall length 115mm
MR-219



Phaco spatula
cum 'Y' rotator
double ended
overall length 145mm
MR-222



Sinsky hook cum ball
manipulator
double ended
overall length 145mm
MR-225



Phaco chop sharp
1mm chopping edge
single ended
overall length 115mm
MR-237



Phaco chop
cum 'Y' rotator
double ended
overall length 145mm
MR-240 with sharp chopper
MR-240A with blunt chopper



Phaco chop blunt
1mm chopping edge
single ended
overall length 115mm
MR-243



Phaco chop
cum spatula
double ended,
overall length 145mm
MR-246 with sharp chopper
MR-246A with blunt chopper



Nucleus spatula for
manipulating, rotating &
cracking the nucleus
during phaco, angled,
overall length 115mm
MR-249



Nucleus paddle tip
spatula for rotating &
cracking the nucleus
during phaco, angled
overall length 115mm
MR-252



Nucleus trisector for
manual small incision
overall length 110mm
MR-255





Divider / Spatula / Irrigating Chopper / 'Y' Rotator / Bimannal

Nucleus bisector for manual small incision surgery, to be used with MR255 overall length 110mm
MR-258



Nucleus cutting board for manual small incision surgery, to be used with MR255 & MR258, overall length 130mm
MR-261



Capsule polishing curette angled overall length 115mm
MR-264



Castroviejo Cyclodialysis Spatula Double ended Ultra Fine overall length 120mm
MR-267



Bi-manual Irrigating hand piece With S.S. 21G 2 irrigating ports overall length 115mm
MR-270



Bi-manual aspiration hand piece S.S. 21G with 0.3mm aspiration port. overall length 115mm
MR-273



Irrigating Nucleus Y-Shaped Manipulator 19G Shaft S.S. dual 0.6mm irrigating port Close to the tip for maximum irrigation. overall length 115mm
MR-276



Irrigating Nuclous Chopper with 0.6mm dual side port 20G overall length 115mm
MR-279 S.S. Gently Curved Shaft
MR-282 S.S. Angled Shaft



Fukasaku Hydro Chopper, blunt distal tip & front notch for engaging iris end. Opening 0.95mm O.D. overall length 115mm
MR-285



Nagahora Hydro Chopper Wedge shaped inferior edge with blunt distal round tip 1.5mm Long ended opening 0.95mm O.D. Overall length 115mm
MR-288



Rosen Hydro Chopper 1.1mm Long, Blunt end 0.95mm O.D. end opening. overall length 115mm
MR-291



Desmarres Retractor
Lid Retractor, solid blade
Small / Med / Large
overall length 140mm
MR-294



Expressor cum
Vectis Double ended
overall length 150mm
MR-297



Expressor Cum
Spatula Double ended
overall length 150mm
MR-300



Lense Loop (Vectis)
overall length 140mm
MR-303



Wilder Lens Loop Serrated
overall length 140mm
MR-306



Catpow Retractor
overall length 130mm
MR-309



Blade breaker swiss made
overall length 115mm
MB-270



Cataract Set with IOL Implantation Set Set of 13 Instruments



1. MF30 Tying forcep straight
2. MF36 Tying forcep curved
3. MF54 McPherson forcep angled
4. MF60 Capsulorhexis Forcep straight
5. MF87 Lims Forcep
6. MF114 Fixation Forcep
7. MF117 Rectus Forcep toothed superior
8. MN312 Barraquer Needle Holder
9. MS240 Corneal Universal Scissor Curved
10. MR204 Repositor Single ended
11. MR207A Repositor cum Spatula double ended
12. MR210 Sinsky Hook extra fine single ended
13. MR219 I.O.L. Manipulator ball ended

Cataract Set Set of 6 Instruments



1. MF 30
2. MF 54
3. MF 87
4. MN 312
5. MS 240
6. MR 210

Cataract Set Set of 4 Instruments



1. MF 30
2. MF 54
3. MF 87
4. MN 312

Phaco Set Set of 4 Instruments





1. MP 30
2. MP 33
3. MR 222
4. MR 243



friMed

Stainless Steel Instruments





Rotator / Phaco Chopper / Pre Chopper

Medelec Phaco spatula cum 'Y' rotator double ended overall size 132mm
MI-03



Medelec Phaco chopper double ended one side sharp & other side blunt
MI-03A



Medelec Phaco chopper, blunt 1.25mm tip length, 1mm chopping edge. overall size 118mm
MI-06



Mrosen Phaco Splitter 60 angled wedge shaped blunt tip.
MI-06A



Medelec Phaco chopper, sharp tip chopping edge 1mm/1.3mm/1.5mm overall size 120mm
MI-09



Medelec Phaco chopper short (0.75mm) pointed "quick chop" tip
MI-09A



Medelec Phaco spatula double ended overall size 142mm
MI-12



Masket combination of the Rosen splitter & a 12mm long round spatula with disc shaped tip
MI-12A



Medelec Phaco chopper cum 'Y' bleichert rotator double ended overall size 130mm
MI-15 with blunt chopper
MI-18 with sharp chopper



Medelec Phaco chopper cum spatula double ended overall size 142mm
MI-21 with blunt chopper
MI-24 with sharp chopper



Pre-Chopper
overall size 95mm
MI-27



Akahoshi Pre-chopper for cracking the soft or hard nucleus Paddle shaped tip which has broad blunt edges on one side to divide a soft nucleus without harming the posterior capsule while the opposite side has the pointed tip and thin edges which are used to penetrate and split a hard nucleus.
MI-27A



Medelec Phaco chopper, angled shaft 'C' shape for entry from the side port, overall size 120mm
MI-30



Medelec Phaco chopper, angled shaft 'S' shape for entry from the tunnel incision overall size 120mm
MI-33



Rosen
Irrigating Nucleus Chopper
with 0.6mm dual side port 20G
MI-33A
MI-33B



Medina Irrigating Nucleus Chopper for micro incision cataract surgery 19G shaft for maximum direct irrigation around the nucleus while chopping sharp / blunt
MI-33C
MI-33D



Irrigating Chopper 19G with notches for iris push for horizontal chopping
MI-33E



Irrigating Sinsky hook for micro incision phaco with dual irrigating ports for manipulating nucleus fragment 20G
MI-33F



Medelec
Irrigating Nucleus Manipulator gently curved 19G shaft dual 0.6mm irrigating port close to the tip for maximum irrigation. The Y-shaped tip for manipulating.
MI-33G



Rosen Wasser Irrigating Endothelial stripper for descemet's Stripping Endothelial keratoplasty
MI-33H





Folders / Inserters

Medelec Folder
for folding Phaco acrylic
soft IOL's through
a small incision
overall size 110mm
MI-36



Medelec Inserter
for inserting Phaco acrylic
soft IOL's through
a small incision
overall size 100mm
MI-39



Medelec Folder
Phaco acrylic soft IOL's
folding forceps with
smooth platform
MI-42



Faulkner
Lens inserting forceps
for inserting acrylic soft IOL
overall size 102mm
MI-45



Medelec Folder
for folding silicon soft IOL
overall size 104mm
MI-48



Medelec Inserter
for soft silicon IOL's inserting
forceps with lock
overall size 115mm
MI-51



Medelec for folding &
holding the shaft IOL's
for insertion
MI-54



Medelec Lens folding forceps
to fold acrylic Soft lens
overall size 130mm
MI-60



Medelec Lens insertion forceps
for small incision to be used
with MI-60 Lens Folder
overall size 110mm
MI-60 A



Medelec Nucleus cracker
overall size 107mm
MI-63



Medelec Nucleus splitter
for cracking the nucleus
overall size 75mm
MI-66



Kansas Nucleus removal
forceps, overall size 109mm
MI-69



Meheta Nucleus cracker
by punching phaco end
non phaco small incision
overall size 150mm
MI-72



Kanas Nucleus Trisector
for manual small incision
overall size 120mm
MI-75



Kanas Nucleus Bisector
for manual small incision
overall size 135mm
MI-78



Kanas Nucleus solid vectis
used as cutting board with
nucleus Trisector and bisector
overall size 138mm
MI-81



Medelec Snare,
nucleus divider
MI-84



Bores Incision spreading
forceps opening jaws
0.5mm
MI-87






Sapphire Ophthalmic Blades

The delicate and sharp brings a tidy edge in a feather touch surgical cut due to high hardness density. Transparency improves microscopic surgical results. The sapphire surgical knife is easily sterilized and retains a strong resistance to corrosion. The cost of per surgery is even lower than using metal knife. Shaughnessy


Microsurgery knife for sideport,
1.0mm straight 45, double bevel
M-10455
2.0mm straight 45, double bevel
M-20455




Phaco slit knife (Keratome)
2.8mm Angled, double bevel
M-2860A
3.0mm Angled, double bevel
M-3060A
3.2mm Angled, double bevel
M-3260A




MVR Blade
1.0mm straight double bevel
MVR-105



Crescent knife
2.8mm Angled single bevel
C-28A
2.0mm Angled single bevel
C-20A



Trapezoid Blade
2.5mm/3.4mm
Four Facet
K-52108



Phaco Knives & Keratomes

Phaco knife for sideport incision
S.S with polycarbonate handle
MI-90



Crescent knife 45 angled 2.8mm
S.S with polycarbonate handle
MI-93



Keratome Angled 45, 2.8mm, 3.0mm, 3.2mm, and
3.5mm. S.S with polycarbonate handle
MI-96



Wound Enlarger St/Angled 45, 4.0mm, 5.1mm,
6.2mm. S.S with polycarbonate handle
MI-99



Speculums

Barraquer
Wire Speculum
MI-102 adult
MI-105 child
MI-108 infant





Barraquer wire speculum extra large
MI-108A



Kratz Barraquer with speculum with open blades for phaco ball ended tips
MI-108B




Kratz-Barraquer Wire Speculum with open blades for unobstructed access to the cataract incision
MI-111 small
MI-114 large



Barraquer Wire Speculum with solid blades
MI-117 small
MI-120 large



Medelec Wire Speculum adjustable spread less pressure
MI-123 soft wire
MI-126 hard wire




Medelec Wire Speculum adjustable spread right / left
MI-126A



Cook Speculum solid blade adjustable spread with lock, also for lasik
MI-129



Speculum Self Retaining
MI-129A




Guyton-Park Eye Speculum with vertical suture posts and canthus hook
MI-132





Speculums

Speculum V-shape
MI-132A



Clark Speculum, clark
MI-135



Speculum adjustable spread
for phaco/lasik overall size 82mm
MI-138



Sauer Lid Speculum
for infants
MI-141



Infant wire speculum
MI-144




Lancaster Speculum, solid blades
MI-147



Temporal Nasal Speculum, adjustable
spread Fenestrated open blades
MI-150




Temporal Nasal Speculum,
adjustable spread solid blades
MI-153




Temporal Approach Wire
Speculum with Spring nasal side
MI-153A
MI-153B
MI-153C



Medelec Barraquer Wire Speculum
with suture notches used to obtain
additional traction for the upper eye lid
MI-153D



Kratz
Aspirating Speculum with open blades
for undistracted access to the cataract
incision with phaco or I/A tip
MI-153E



Lieberman
Aspirating speculum, adjustable for
unobstructed access to the cataract
incision with Phaco or I/A tip
MI-153F



Scleral Depressors / Retractors

Scleral Depressor
with thimble
MI-156



Scleral Depressor double ended
overall size 130/128mm
MI-159



Localizer & Scleral
depressor for fixation & rotation
overall size 127mm
MI-162



Sheet Glide Intraocular lens glide
MI-165



Jaffe wire lid retractor may be clamped
to surgical drape with suture or elastic band
overall size 57mm
MI-168



Stallard ptosis plate
overall size 52mm
MI-171



Desmarres
Lid Retractor, solid blade
size-1 - overall size 135mm
size-2 - overall size 140mm
size-3 - overall size 148mm
size-4 - overall size 150mm
MI-174



Jaeger Lid plate
overall size 100mm
MI-177



Tenant Lid retractor
overall size 24mm
MI-180



Cat Paw Lacrimal
retractor 4-prong
overall size 135mm
MI-183





Retractors / Calipers Fixation Ring

Mueller Lacrimal
Sac retractor
MI-186



Stevenson Lacrimal Sac retractor
3x3 prong, self retaining
MI-189



Dastoor Iris retractor
straight / right / left
MI-192



Schepens Arruga Retractor
overall size 155mm
MI-195



Kilner Retractor
for Occuloplasty
MI-198



Thorton Globe fixation ring
multiple teeth 16mm dia
MI-201



Fixation Ring for phaco &
IOL insertion for clear cornea
overall size 115mm
MI-204



Globe fixation ring with
concentric grooves for
stabilizing the globe during
PRK & Lasik procedures.
MI-204A



Globe fixation ring with concentric
grooves for stabilizing the globe &
providing counter pressure for clear
cornea phaco incisions & IOL insertion
MI-204B



Fleringa Scleral fixation ring
size 12mm to 22mm
MI-207




Scleral ring for global
support during keratoplasty
set of three
MI-207A



Castroviejo Caliper, 20mm spread
single/double graduated
overall size 75mm
MI-210 straight
MI-213 curved



Stahl A/C Caliper for
measuring corneal diameter
overall size 36mm
MI-216



Rule Stainless steel
overall size 150mm
MI-219



Trephines

Castroviejo Corneal Trephine
 Size in mm 5, 5.5, 6, 6.5, 7, 7.5,
 8, 8.5, 9, 9.5, 10.0, 10.5, 11.0, 11.5
 & 12mm overall size 74mm
MI-222



Disposable Corneal Trephine
 Size in mm 5, 5.5, 6, 6.5, 7, 7.5,
 8, 8.5, 9, 9.5, 10.0, 10.5, 11.0, 11.5
 & 12mm overall size 42mm
MI-225



Elliot Scleral Trephines
 Size in mm 1.0, 1.5 & 2.0mm dia
MI-228



Arruga bone Trephines with
 sharp cutting edge with teeth,
 used for cutting the bone
 Size in mm 7, 7.5, 8, 8.5, 9 & 10mm
MI-231 without center drill
MI-234 with center drill



R. K. Instruments

Corneal / Radial Marker
 4 / 6 / 8 / 12 / 16 Blades
 overall size 105mm
MI-237



Mendoz Degree Gauge for
 accurate axial alignment of
 corneal incisions to correct
 12mm I.D.
MI-240



Bores Axis Marker used with
 mendoz Degree gauge to
 mark axis in which corneal
 incision for correction of
 astigmatism are to be made
MI-243



Bores T-Incision Marker has two
 6mm long marking blades with spread
 overall size 103mm
MI-246



Hoffer Optic zone marker with
 cross hair 3.0mm to 9.0mm
 (in variation of 0.25mm)
 overall size 102mm
MI-249



Lasik optical zone marker
 with spatula for flap size
 6.5mm, 7mm, 7.5mm & 8mm
 overall size 115mm
MI-252





Curettes Polishers

Meyerhoefer
Chalazion curette
size 1 / 2 / 3 & 4
overall size 130mm
MI-255



Nightingale capsule polishing
curette 1.75 mm dia open ring
with semi sharp edge
overall size 120mm
MI-258 posterior
MI-261 anterior



Anis Capsule Formix Scrubber
gently curved flexible shank with
sand blasted ball shaped tip
overall size 112mm
MI-264




IOL Hooks

Jaffe-Knolle
Iris Hook, small angled
overall size 125mm
MI-267 angled
MI-270 straight



Sinsky Lens
manipulating hook blunt tip
overall size 118mm
MI-273 angled
MI-276 straight



Medelec
Sinsky Hook
Reversed hook angled
with 0.2mm blunt tip
MI-276A



Kuglen Iris hook and
lens manipulator
push-pull model
MI-279 angled
MI-282 straight



Hirschman
Lens spatula 0.5mm
wide with small notches
for manipulating IOL's
MI-285



Manipulators

Jaffe-Beachert
Lens Nucleus
rotator, angled
overall size 128mm
MI-288



Clayman Nucleus
rotator V-notched tip
MI-291



Lester no hole
IOL manipulator
overall size 118mm
MI-294



Drysdale Nucleus
manipulator paddle
shaped tip
MI-297



AKAHOSHI Nucleur
Sustainer Micro ring
tip angled shaft to be
used with prechopper
MI-297A



Medelec Sinsky hook cum
push-pull ball manipulator
double ended
overall size 147mm
MI-300



Medelec Sinsky hook cum
micro spatula double ended
overall size 147mm
MI-303



Medelec Sinsky hook cum
chopper blunt / sharp double
ended, overall size 136mm
MI-306



Hunkeler hemi-ball angled shaft
for dialing lenses without holes
MI-309



Fenzl lens manipulator
"T" shaped, angled with
0.2mm blunt tip
overall size 120mm
MI-312



Medelec
End Capsular Controller
trumpet shaped tip
MI-312A





Spuds / Loops / Cystotomes

Medelec Foreign body spud
MI-315



Medelec Needle & Spud in sliding case
overall size 165mm
MI-318



Lewis Lens loop (vectis)
small / large
overall size 142mm
MI-321



Wider Lens loop serrated
overall size 138mm/142mm
MI-324



Simcoe Lens nucleus loop,
fine serrations, strong curve
double ended
MI-327



Cystotome and spoon
MI-330



Spoons / Spatulas / Repositors

Bunge
evisceration spoon
MI-333 small 8mm
MI-336 large 12mm



Mule evisceration spoon
overall size 140mm
MI-339



Wells enucleation spoon
overall size 140mm
MI-342



Paton spatula & spoon
overall size 140mm
MI-345



Cullers iris spatula
sterling silver blade
overall size 122mm
MI-348



Barraquer cyclodialysis spatula
single ended 1mm wide
overall size 120mm
MI-351



Castroviejo cyclodialysis spatula
double ended, delicate
overall size 140mm
MI-354



Spoons / Spatulas / Repositors

Castroviejo cyclodialysis spatula flat double ended overall size 140mm
MI-357



Lens spatula double ended overall size 160mm
MI-360



Medelec Iris repositor double ended overall size 145mm
MI-363



Medelec Repositor curved / angled double ended, overall size 142mm
MI-366



Medelec Nucleus spatula for manipulating, rotating & cracking the nucleus during phaco
MI-369



Membrane Peeler Set



Membrane Peeler
MI-372



Membrane Spatula curved
MI-372B



Membrane Spatula hooked with pointed tip
MI-375A



Membrane Scratcher claw type
MI-378



Membrane Spatula hooked
MI-372A



Membrane Spatula hooked short length
MI-375



Membrane Scratcher hooked
MI-375B



Membrane Spatula knob
MI-378A

Knives

Barraquer Cataract knife in sliding case
MI-381



Ziegler's Knife needle 1.5 cutting edge overall size 132mm
MI-384 in sliding case
MI-387 without sliding case



Paufique Graft knife angled tip, curved cutting edge for keratoplasty overall size 108mm
MI-390



Tooks Corneal knife, straight blade with curved cutting edge overall size 136mm
MI-393




Gill Corneal knife, curved cutting edge
MI-396



Barraquer keratoplasty knife sharp pear shaped for lamellar dissection
MI-399



Hockey stick shaped knife for peeling & folding back the epithelial flap
MI-399A





Muscle Hooks / Cauteries / Expressors

Graefe strabismus hook
size 1,2,3
overall size 150mm
MI-402



Jameson muscle hook
small / large
overall size 130mm
MI-405



Smith lens expressor
overall size 150mm
MI-408



Kirby expressor cum spatula double ended
overall size 147mm
MI-411



Kirby hook cum vectis double ended
overall size 145mm
MI-414



Kirby expressor cum vectis double ended
overall size 155mm
MI-414A



Stevens tenotomy hook
overall size 132mm
MI-417



Chavasse strabismus hook curved shaft
MI-420




Culler Rectus muscle hook
MI-423



Gass retina detachment hook
MI-426



Wills cautery
s.s.with copper
MI-429 ball point size 135mm
MI-432 flat point
MI-435 st point size 140mm



Blade Breakers / Dissectors

Swiss model blade breaker
and holder small / big
overall size 95mm/98mm
MI-438



Barraquer blade breaker and holder
overall size 102mm
MI-441



Blade cutter , ½ size
MI-444



Dastoor lacrimal sac dissector
double ended with sharp
gauge and blunt dissector
overall size 132mm
MI-447



Lang lacrimal sac dissector
double ended with sharp
dissector and curette
overall size 140mm
MI-450



Bone Punches / Nibblers

Kerrison
bone punch 85mm sheft
size 1.5mm, 2.0mm, 3.0mm & 4.0mm
MI-453



Citeli
bone nibbling rongeur
80mm sheft
MI-456



Hajeks
bone nibbling rongeur 105mm sheft
MI-459 upward cutting
MI-462 downward cutting





Bone Punches / Nibblers

Beyur
bone rongeur
overall size 165mm
MI-465 single action
MI-468 double action



Kelly (Glaucoma) Punch
overall size 140mm
MI-471




Holth
corneoscleral punch
1.5mm & 2mm bite
MI-474



Lacrimal Instruments

West Lacrimal chisel,
straight, overall size 150mm
MI-477



West Bone gauge
overall size 150mm
MI-480



Freer Periosteal
elevator double
ended, overall size 152mm
MI-483



Traquair Periosteal
elevator double ended
with blunt ended tips
overall size 150mm
MI-486



Mallet
overall size 175mm
MI-489



Nasal speculum
Adult / Child
MI-492



Tilly nasal packing forcep
MI-495



Probes / Dilators

Harms Trabeculotomy probe
right / left ,overall size 50mm
MI-498




Pigtail probe with suture holes
overall size 155mm
MI-501



Wilder Lacrimal dilator
overall size 110mm
MI-504 No:1, fine taper
MI-507 No:2, medium taper
MI-510 No:3, heavy taper



Bowman Lacrimal probe,
malleable stainless steel,
overall size 140mm/145mm/150mm
No:- 0000-000, 000-00,
00-0, 1-2, 3-4 and 5-6
MI-513



Catalano Lacrimal set
of two blunt s.steel malleable
probes with silicon tubing
stainless steel
MI-516
MI-516A Incubation 23G



SQ Handle Forceps / Scissors

Medelec
Squeeze handle scissor i
ntraocular & vitreous in 20G
Overall size 150mm
MI-519 vertical 90
MI-522 horizontal opening straight
MI-522A horizontal curved
MI-522B horizontal angled 45



Medelec Squeeze handle forceps
intraocular & vitreous in 20G
overall size 150mm
MI-525 vertical
MI-528 horizontal
opening straight



Wilson foreign body forceps
three 2mm prongs in
20 gauge shaft
MI-531





Clamps / Forceps

Serrefine Bulldog clamp
small straight / curved
overall size 55mm
MI-534



Baby jones Towel clamp
cross action overall size 75mm
MI-537



Backhouse Towel clamp
overall size 85mm
MI-540



Hartman Mosquito forceps
delicate, 33/4" (95mm) long
straight and curved
MI-543



Halsted Mosquito forceps
delicate, 53/4" (145mm) long
straight and curved
MI-546



Botvin Iris forceps 1x2 teeth
for peripheral iridectomies
overall size 73mm
MI-549



Hess Iris forceps 2mm teeth
standard model for peripheral
iridectomies
overall size 75mm
MI-552



Kirby Iris forceps 2mm teeth
MI-555



Colibri forceps, delicate 1x2 teeth
0.4mm with tying platform
overall size 75mm
MI-558



Polack Corneal Forcep,
very delicate 1mm spread
1x2 teeth 0.12mm
MI-558A



Hoskin colibri
corneal forcep very delicate
with platform 1x1 pierse 0.1mm
overall size 75mm
MI-561



Medelec Hoskin corneal forcep
castroviejo long wide handle
with platform, overall size 100mm
overall size 100mm
MI-561A




Forceps

Pierse corneal forceps very delicate pierse teeth 0.1mm with tying platform overall size 88mm
MI-564 straight
MI-567 curved



Bonn Iris forceps 1x2 teeth overall size 80mm
MI-570



Mc. Pherson corneal forceps 1x2 teeth extra fine with tying platform
MI-573 long overall size 100mm
MI-576 medium overall size 100mm



Medelec Corneal forceps straight micro 1x2 teeth 0.12/0.2mm with tying platform overall size 88mm
MI-579



Kirby Corneal forceps 1x2 teeth, fine, with tying platform
MI-582 straight, size 98mm
MI-585 curved, size 95mm



Girard Corneoscleral forceps 1x2 teeth 0.12mm with tying platform overall size 104mm
MI-588




Kelman-McPherson forceps 8mm long angled jaw 1x2 teeth 0.3mm overall size 85mm
MI-591



Lims 1x2 teeth extra fine, forcep with tying platform overall size 92mm
MI-594



Medelec Lims 1x2 extra fine wide handle forcep with tying platform overall size 90mm
MI-597



Castroviejo corneal wide handle forceps 1x2 teeth extra fine with tying platform, overall size 110mm
MI-600 straight
MI-600A curved



Bores corneal fixation forceps U-shaped 3mm spread 1x2 teeth at each tip 0.12mm
MI-603 straight
MI-606 angled





Forceps

Kramer corneal fixation forceps with lock U-shaped 8mm or 13mm spread straight / angled overall size 110mm
MI-609



Jayles fixation cum corneal 'U' forceps 1x2 teeth with 6mm tying platform overall size 110mm
MI-612



Kelman-McPherson forceps angled very delicate 7.5mm long smooth jaws overall size 105mm
MI-615



Kelman-McPherson forceps, angled very delicate
MI-618 Short model, 82mm
MI-618A long model, 110mm



Kelman-McPherson lens holding forceps, angled very delicate 11mm long smooth jaws overall size 105mm
MI-621



Medelec IOL forceps angled very delicate, 10mm long smooth jaws
MI-624



Kraff capsule tag forceps, 11mm long jaws with 1.5mm long smooth platform at tips overall size 105mm
MI-627



Jaffe tying forceps very delicate, smooth jaws with 6mm long platform for 10.0 & 11.0 sutures
MI-630 straight overall size 100mm
MI-633 curved overall size 102mm




Castroviejo tying forceps, delicate with platform, straight overall size 83mm
MI-636 straight
MI-639 curved




Forceps

Castroviejo tying forceps, with platform serrated wide handle overall size 110mm
MI-639A straight
MI-639B curved




Barraquer tying forceps with 7mm tying platform for 8.0mm sutures overall size 112mm
MI-642 straight
MI-645 curved




Mcpheerson tying forceps 5.0mm long tying platform overall size 85mm
MI-648 angled



Tenant tying forceps with guide pin extra delicate with 6mm long platform for 9.0 to 11.0 sutures overall size 100mm
MI-654 straight
MI-657 curved



Anis tying forceps extra delicate smooth jaws for 9.0 to 11.0 sutures
MI-660



Sinsky tying forceps, very delicate smooth jaws for 8.0 to 11.0 sutures
MI-663



Medelec Capsulorhexis forceps Ultrafine, new model overall size 86mm
MI-666



Medelec Capsulorhexis forceps very delicate overall size 104mm
MI-669



Medelec Capsulorhexis Forcep 23 gauge Curved Shaft with micro grasping tip
MI-672



Capsulorhexis Forcep 21/23G Curved Shaft
MI-672A





Forceps

Utrata Capsulorhexis forceps
shanks ,straight / curved
overall size 94mm
MI-675



Masket Capsulorhexis
forceps very delicate tip,
gently curved
overall size 104mm
MI-678



Inamura Caprulothexis
Forcep delicate
MI-678A



Dressing forceps, serrated
overall size 100mm
MI-681 straight
MI-681A curved



Bonaccolto conjunctival forceps
longitudinally serrated jaws
with cross serrations at tips
overall size 100mm
MI-684



Moorfield conjunctival forceps
overall size 100mm
MI-687




Tissue fixation, 1x2 teeth
small / large
overall size 100mm
MI-690



Stern Castroviejo Fixation forcep
1x2 teeth with lock straight
MI-690A



Moody Fixation forcep
1x2 teeth with lock
MI-690B - Left
MI-690C - Right



Fixation forceps
1x2, 2x3, and 3x4 teeth
overall size 110mm
MI-693



Graefe fixation forceps
MI-696 7/8 teeth, 100mm
MI-699 2x3/3x4 teeth, 115mm



Green fixation forceps
multiple teeth
overall size 100mm
MI-702 with lock
MI-705 without lock



Dastoor superior rectus
forceps 1x2 teeth
overall size 115mm/105mm
MI-708



Forceps

St. Martin fixation cum corneal forceps 1x2 teeth with tying platform overall size 88mm
MI-711



Cilea forcep
MI-711A



Beer cilea forceps overall size 90mm
MI-714



Barraquer cilea forceps smooth jaws straight / slightly angled
MI-717



Jewelers forceps extra delicate overall size 110mm
MI-720 straight
MI-723 curved



Jewelers forceps extra fine overall size 105mm
MI-726



Arruga capsular forceps standard curve overall size 102mm
MI-729



Arruga capsular forceps new curve overall size 100mm
MI-732



Arruga capsular forceps cross action
MI-732A



Shepard lens holding forceps with serrated lower jaw gently curved shanks overall size 120mm
MI-735



Kratz lens holding forceps very delicate, curved narrow smooth jaws overall size 112mm
MI-738



Clayman lens holding forceps delicate angled jaws gently curved shanks overall size 122mm
MI-741





Forceps

Sinsky lens holding forceps,
delicate T-shaped lower jaw,
suitable for ridged lenses
overall size 120mm
MI-744



Moore lens holding forceps
overall size 118mm
MI-747



Blayds lens holding forceps
overall size 110mm
MI-750



Medelec lens holding forceps,
McPherson type overall size 112mm
MI-753 long blade
MI-756 short blade



Daljit lens holding forceps
overall size 90mm
MI-759



Jamson muscle forceps
with side lock with teeth,
adult size right / left
overall size 100mm
MI-762



Worth muscle forceps
11mm longitudinal serrated
jaws with sliding lock
overall size 95mm
MI-765 right
MI-768 left



Medelec muscle cum suture
holding forcep cross action
MI-768A



Berke ptosis clamp
with side lock
overall size 100mm
MI-771



Snellen entropion forceps
right / left overall size 98mm
MI-774



Desmarres chalazion forceps, oval large
20mm, 26mm, 31mm overall size 95mm
MI-777



Hunt chalazion forceps round 12mm diameter
overall size 95mm
MI-780



Forceps

Lambert chalazion forceps round
10mm/15mm diameter, overall size 95mm
MI-783



Ayer chalazion forceps
overall size 95mm
MI-786



Alabama Utility Forcep straight shaft
with 15mm lone 'V' grooved platform
MI-789A



Fechtner Conjunctival Suture
forceps with delicate rings to
hold the conjunctiva without
tearing with tying platform
overall size 112mm
MI-789B



Scissors

Castroviejo corneoscleral scissors
small / medium blades
7mm/10mm blades
overall size 95mm
MI-792 right
MI-795 left
MI-798 R/L with lock



Castroviejo corneal scissors, micro
universal curved, fine blunt tips
overall size 110mm
MI-801 small blades
MI-804 medium blades
MI-807 large blades



Castroviejo micro corneal scissors
universal curved serrated body,
11mm blades, overall size 95mm
MI-810 small blades





Scissors

Medelec corneal transplant scissors, strong curve overall size 120mm
MI-813



Castroviejo micro corneal scissors angled on flat, 11mm blades overall size 108mm
MI-816 small blades



Iris scissors, sharp tips small blade, straight / curve overall size 115mm
MI-819



Fine stitch scissors sharp tips, medium blade overall size 120mm
MI-822



Westcott tenotomy scissors curved blunt tips overall size 122mm
MI-825



Westcott Stitch Scissor very sharp pointed tips, standard blades, curved
MI-825A



Westcott tenotomy scissors curved sharp pointed tips overall size 120mm
MI-828



McPherson- Westcott Stitch Scissor very sharp pointed tips, small blades, curved
MI-828A



Castroviejo keratoplasty scissors angled to side, blunt tips overall size 110mm
MI-831 small blade
MI-834 long blade



Scissors

Vannas straight scissors,
sharp tips,
overall size 90mm
MI-837



Vannas curved scissors
sharp / blunt tips,
overall size 93mm
MI-840



Medelec Vannas scissor Kirby
round handle curved blade
overall size 95mm
MI-840A



Gills-Vannas capsulotomy scissors,
extra fine long blades, angled upward
overall size 95mm
MI-843



Medelec Capsulotomy scissor Kirby
round handle angled upward
overall size 87mm
MI-843A



Vannas scissors, angled to side
overall size 80mm
MI-846



Medelec capsule tag cutting
scissors, curved blunt round tips
overall size 100mm
MI-849



Stern-Gills scissors, extra fine
angled upward / straight
overall size 105mm
MI-852




Barraquer Iris scissors, blunt tips
7mm/10mm blades
overall size 57mm
MI-855



De-weckers Iris scissors
12mm blades, sh/sh, bl-bl
overall size 110mm
MI-858




Eye scissors delicate pointed
iris scissors
overall size 90mm
straight / curved
MI-861





Scissors

Eye scissors
conjunctival sharp, pointed tips
105mm straight / curved
MI-864



Stevens
tenotomy scissors, standard
blunt tips, straight / curved
MI-867



Medelec strabismus scissors
rounded blunt tip
straight / curved
overall size 100mm
MI-870



Medelec stitch scissors delicate
straight / curve
overall size 100mm
MI-873



Enucleation scissors
MI-876 slight curve
MI-879 medium curve
MI-882 full curve



Needle Holders

Barraquer needle holders,
very delicate jaws straight
overall size 130mm
MI-885 with lock
MI-888 without lock



Barraquer needle holders,
very delicate jaws curved
overall size 115mm
MI-891 with lock
MI-894 without lock



Barraquer needle holder, micro
curved 10mm long delicate jaws
overall size 115mm
MI-897 without lock
MI-900 with lock



Needle Holders

Barraquer Troutman
needle holder, micro
strongly curved 10mm
long jaws without lock
overall size 112mm
MI-903



Medelec Needle holder
Barraquer long, micro delicate
overall size 110mm
MI-906
MI-906A straight
MI-906B curved



Castroviejo micro needle holders
curved small 3 1/2"
overall size 90mm
MI-909



Castroviejo micro needle
holders curved medium
overall size 118mm
MI-912



Castroviejo needle
holders, long straight
overall size 140mm
MI-915 with lock
MI-918 without lock



Castroviejo needle holders,
long curved
MI-921 with lock, size 130mm
MI-924 without lock, size 140mm





Needle Holders

McPherson needle holders
 very delicate jaws
 overall size 95mm
MI-927 straight
MI-930 curved



Kalt needle holder for
 lids and superior rectus
 suture with locking handle
 overall size 135mm
MI-933



Arruga needle holder
 with locking handle
 overall size 138mm
MI-936



Silcock needle holder for
 lid and superior rectus
 suture with locking handle
 overall size 140mm
MI-939



Cannulas

Air Injection cannula, angled 5mm from tip
MI-942



Healon aspirating cannula
 angled 10mm from tip
MI-945




Rainin air injection cannula
 spatulated tip, angled at 7mm
MI-948



A.C. Wash Cannula
 18,20,21 & 22 gauge
MI-951



Mc. Intyre anterior chamber irrigating
 cannula, smooth blunt tip, 20 to 27G
MI-954



Randolph cyclodialysis cannula
 curved spatulated tip, angled at 13mm
MI-957 front opening



Lasik cannula 25g angled tip
 with 6 irrigating parts
MI-960



Cobra
 Lasik Cannula/Vectis 6 Ports
 for multi directional irrigation 23G
MI-960A



Cannulas

Lacrymal Cannula
half/ full curved 23-26G
MI-963



Welsh flat olive tip cannula 25G
straight / left / right
MI-966



Kelman cystotome with sharp
point and sharp cutting edge 22G
straight / angled
MI-969



Jenson posterior capsule polisher
all round sand blasted bulbous tip
MI-972



Greather collar button, micro iris
retractor & posterior capsule polisher
with side irrigation straight / angled
MI-975



Simcoe cortex extractor cannula,
smooth closed tip, with top / side
opening facing into curve 23G
MI-978



Knolle-Pearce irrigating vectis dual
fluid flow through three front
opening ports 23G loop
MI-981 short
MI-984 long



Nishi Nuclis Delivery
Cannula with irrigating port
MI-987



Alfonso-Mc.Intyre
Irrigating Nucleus spoon single 0.3mm
front opening port for extracting nucleus
fragment during manual small incision surgery
MI-987A



Simcoe Nuclis Delivery Cannula right-left
two front opening ports, serrated
MI-990




Sheet irrigating vectis dual fluid flow
through three large front opening
ports, serrated 20 gauge loop
MI-993



Nucleus Delivery Punch
with tapered plate
MI-996



Medelec
Irrigating Vectis Serrated
with cutting punch in center
MI-996A





Cannulas

Halfman
Nucleus loop, 3mm wide
MI-999



Medelec
Irrigating Vectis narrow fine long
MI-999A serrated




Sinsky lens manipulating
cannula front opening, angled
MI-1002



Agarwal Irrigating Nucleus chopper
semi-sharp dual side ports 20G
MI-1002A



A.C. Maintainer
with silicon tubing
threaded on tip
MI-1005



Medelec Infusion cannula
20G beveled tip with luer lock
adopter for vitreo retinal surgery
MI-1005A



Nucleus Hydrodissecting
cannula, 24G, 25G, 26G
MI-1008



Hydrolysis cannula
MI-1011



Medelec Hydro dissection Cannula
flat tip angled 90 with beveled
opening 26/27G
MI-1011A



J-shaped cannula
MI-1014



J-Loup cannula,
irrigating / aspirating
with silicon tubing
right / left
MI-1017





Cannulas

Simcoe (Regular) Direct

Irrigating-Aspirating cannula with silicon tubing. Original model, with side aspiration through tubing hub and irrigating through luer lock hub. 23gauge with 0.3mm dia aspiration port facing into curve and slightly shorter irrigation tip with 45 angled, front opening. Die punch hub

MI-1020

MI-1023 20G, 21G & 22 Gauge available





Reversed Simcoe

I/A cannula, side aspiration through luer lock hub, shorter end irrigation through tubing hub.

Available : 20G, 21G, 22G & 23G Dia Punch Hub

MI-1026



Simcoe Irrigation & aspiration Cannula with serrated shaft for Hydro expression of nucleus for non-phaco small incision surgery

MI-1026A



Gills Irrigating Aspirating Cannula front opening with irrigating port 0.5mm from tip angled, 23G

MI-1026B



Silicon Bulb with adopter
MI-1029



Daljit
I/A Cannula with
Silicon tubing
MI-1032





Cannulas

Infusion handle
small
MI-1035



Irrigating-Aspirating Handle
stainless steel, long
MI-1038



Bard-Parker Handle
MI-1041




Hand Piece with silicon
tubing self locking with
cortex aspirating cannula
MI-1041A



Aspirating Hand piece with 21G
0.3mm aspirating port and fitted
with silicon tubing
MI-1041B



Medelec
Bimanual Irrigating
Hand Piece 21G with
two irrigating port
MI-1044



Medelec
Bimanual Aspiration
Hand Piece 21G with
0.3mm aspiration port
MI-1047



Medelec Instrument Caddy
made of stainless steel for sterilization
of micro surgical instruments, with
MI-1050 single soft silicon matt. Size 8"x4 1/2
MI-1053 with double silicon matt, big size



Sterilizing case autoclavable plastic with silicon matt
MI-1053A medium
MI-1053B large
MI-1053C double level includes removable tray & 2 silicon matts



MI-1053A



MI-1053B



MI-1053C

LASIK Flap
Protector & Spatula
MI-1500



LASIK Flap Protector blades
angled 45° for retracting & protecting
superiorly hinged corneal flaps.
MI-1503



LASIK Flap Protector blades
angled 15° for retracting and
protecting superiorly
hinged corneal flaps.
MI-1506



LASIK Spatula dual 12mm
long 0.25mm diameter arms
MI-1509



LASIK Flap Edge
Dissector For LASIK
enhancements.
MI-1512



LASIK Flap Manipulator
delicate olive tip
MI-1515



Stabilization Ring with swivel
handle, 13mm dia. with blunt
teeth on both sides For
stabilizing either the right
or left eye during LASIK.
MI-1518



LASEK Trephine
270° semi-sharp edge,
90° cutout for flap hinge
MI-1521 (8/9mm diameter)



LASEK Trephine &
Alcohol Retaining Well
MI-1524 (8mm trephine with 8.5mm well)
MI-1524A (9mm trephine with 9.5mm well)





Lasik / Lasek

LASEK Alcohol Well for retaining alcohol solution on the corneal epithelium.

MI-1527 (8.5mm diameter)

MI-1527A (9.5mm diameter)



LASEK Epi-Dissector spatulated tip with semi-sharp edges for undermining the epithelium.

MI-1530



LASEK Micro Hoe semi-circular tip with beveled edge For lifting the precut margin of the epithelial flap.

MI-1533



LASEK Epi Peeler semi-sharp edges for peeling and folding back the epithelial flap.

MI-1536



LASEK Micro Hoe & Epi Peeler

MI-1539



Corneal Marker for LASIK surgery for identification & realignment of the lamellar flap.

MI-1542



Corneal Marker for LASIK surgery For identification & realignment of the lamellar flap.

MI-1545



LASIK spatula & flap unzipper overall size 144mm

MI-1548



LASIK spatula double ended overall size 143mm

MI-1551




Lasik spatula single ended

MI-1551A



Aspirating speculum, adjustable for unobstructed access to the cataract incision with Phaco or I/A tip
MI-1554



Speculam adjustable spread for phaco/lasik, overall size 82mm
MI-1557

Multi-Purpose LASIK Forceps spatulated jaws & tips smooth grasping surfaces, vaulted for lifting, dissecting & grasping the corneal flap.
MI-1560



Buratto lasik flap holding forceps with fine disk & serration for firm grip
MI-1560A



Spatulated LASIK Cannula flattened spatulated tip end opening
 26 gauge
MI-1563



LASIK Cannula with one front and two side ports
MI-1566



LASIK Cannula low profile tip end opening
 23 gauge
MI-1569



LASIK Fountain Cannula smooth bullet shaped tip single centre port for directing fluid toward the superior hinge
MI-1572



LASIK Cannula bullet shaped tip 6 ports for multi-directional irrigation
MI-1575



LASIK Cannula closed spatulated tip 6 ports for multi-directional irrigation
MI-1578



Buratto LASIK Cannula one front & two side ports 23G
MI-1581



LASIK Cannula double-armed, 8 irrigating ports
MI-1584




LASEK Alcohol Cannula delicate olive shaped tip
 25G
MI-1587







① ⑤ ⑥ ⑦ ⑧ ⑨ ⑩ ⑪ ⑫ ⑬




⑭ ⑮ ⑯ ⑰ ⑱ ⑲ ⑳ ㉑ ㉒



㉓ ㉔ ㉕ ㉖ ㉗ ㉘ ㉙ ㉚ ㉛ ㉜ ㉝






㉞ ㉟ ㊱ ㊲ ㊳ ㊴ ㊵ ㊶ ㊷ ㊸ ㊹



㉚ ㉛ ㉜ ㉝ ㉞ ㉟ ㊱ ㊲ ㊳ ㊴

MS1059. Intra/extra Capsular Cataract With IOL Implantation Set

1. MI03 Spatula cum Y rotator double ended
2. MI102 Barraquer wire speculum adult
3. MI105 Barraquer wire speculum child
4. MI108 Barraquer wire speculum infant
5. MI192 Iris wire retractor straight
6. MI210 Caliper Graduated
7. MI258 Capsule polishing Curette
8. MI273 Sinsky hook angled
9. MI279 Kuglen hook and manipulator angled
10. MI321 Vectis/loop
11. MI357 Cyclodialysis spatula double ended
12. MI360 Lens spatula double ended
13. MI363 Iris repositor double ended
14. MI381 Cataract knife sliding mode
15. MI393 Corneal tools knife
16. MI408 Smith Lens expressor
17. MI429 Ball heat Cautery
18. MI438 Blade breaker and holder
19. MI534 Serrafine bulldog clamp 2nos
20. MI537 Baby jones towel clip cross action
21. MI543 Hartman mosquito forceps straight
22. MI546 Halstead mosquito forceps curved
23. MI549 Botvin Iris forceps
24. MI564 Pierson-Hoskin corneal forceps straight
25. MI579 Corneal forceps toothed straight
26. MI594 Lims toothed forceps extra fine
27. MI621 Kelman McPherson IOL forceps
28. MI630 Tying forceps straight
29. MI633 Tying forceps curved
30. MI681 Dressing forceps serrated jaws
31. MI687 Moorfield Conjunctival forceps
32. MI690 Fixation forceps
33. MI708 Superior Rectus forceps
34. MI729 Capsular Arruga forceps
35. MI735 Lens holding forceps
36. MI804 Castroviejo corneal scissor universal
37. MI810 Medelec Micro Corneal Scissor
38. MI825 Westcott tenotomy scissors
39. MI834 Corneal Scissors angled to side
40. MI837 Vannas Micro Scissors straight
41. MI843 Vannas Capsulotomy Scissor angled Upward
42. MI894 Barraquer Micro Needle Holder
43. MI912 Castroviejo Needle holder standard
44. MI942 Air injection Cannula
45. MI972 Posterior capsule polisher
46. MI978 Cortex cannula
47. MI1020 I/A Simcoe Cannula
48. MI1029 Silicon bulb with connector
49. MI1038 Infusion handle
50. MI1041 B.P.Handle
51. MI1050 Instrument caddy with silicon matt




Set-2 Ms1062. Non-Phaco Set (Small Incision)

1. MI63 Nucleus Cracker
2. MI69 Nucleus Removal forceps
3. MI75 Nucleus Trisector for manual Non-phaco
4. MI78 Nucleus Bissector
5. MI81 Nucleus cutting board
6. MI90 Side port knife disposable
7. MI93 Crescent blade disposable
8. MI96 Keratome blade disposable
9. MI123 Speculum self retaining
10. MI369 Nucleus spatula for cracking, rotating
11. MI438 Blade breaker and holder
12. MI588 Girard Corneoscleral forceps
13. MI609 Kramer fixation U-shaped 1x2 teeth.
14. MI744 Lens holding forceps T-shape
15. MI969 Bluementhol Cystotome
16. MI978 Cortex Cannula
17. MI993 Irrigating vectis serrated
18. MI996 Nucleus delivery punch cannula
19. MI1005 A.C.Maintainer with tubing
20. MI1008 Nucleus Hydrodissection cannula
21. MI1011 Hydrolysis cannula angled
22. MI1017 J-loup cannula I/A RT/LT-2nos
23. MI1047 Medelec asperation hand piece
24. MI1050 Instrument Caddy with silicon matt

Set-3 Ms1062A. Non-Phaco Set (Small Incision) with Sapphire surgical Blade

25. Sapphire Surgical Blade Set Of Three
 - A. M1045S
 - B. M3060A
 - C. C28 A



Phaco Emulsification Set



1. MI03 Phaco spatula cum Y rotator double ended
2. MI06 Phaco chop blunt
3. MI09 Phaco chop sharp
4. MI27 Pre-choper
5. MI36 Lens folder acrylic soft
6. MI39 Lens inserter acrylic soft
7. MI48 Medelec folder for silicon soft lens
8. MI51 Medelec inserter for silicon soft lens
9. MI63 Nucleus cracker
10. MI69 Nucleus removal forceps
11. MI138 Phaco speculum self retaining
12. MI153 Temporal nasal speculum
13. MI204 Phaco fixation ring
14. MI210 Castroviejo caliper graduated
15. MI297 Nucleus manipulator paddle shape
16. MI354 Cyclosidalysis spatula double ended
17. MI369 Nucleus spatula for cracking, rotating
18. MI534 Serafine Bulldog clamp (2)
19. MI537 Towel clip (2)
20. MI561 Hoskin Colibri forceps
21. MI594 Lims forceps toothed extra fine
22. MI609 Kramer corneal fixation forceps
23. MI630 Tying forceps straight
24. MI666 Capsulorhexis Ultrafine forceps
25. MI843 Gills vannas capsulotomy scissors
26. MI897 Needle Holder micro
27. MI942 Air Injection cannula
28. MI1008 Hydrodissection cannula
29. MI1011 Hydrolysis cannula
30. MI1029 Silicon bulb with connector
31. MI1044 Bimanual irrigation hand piece
32. MI1047 Bimanual aspiration hand piece
33. MI1053 Sterilizing case with matts



Lacrimal / D.C.R. Set



1. MI129 Speculum cook adjustable solid blade
2. MI174 Lid retractor medium
3. MI183 Cat Paw retractor (2)
4. MI186 Mueller Retractor
5. MI234 Bone drill with/without centre drill (2)
6. MI438 Blade breaker and holder
7. MI447 Lacrimal sac dissector double ended
8. MI453 Karrison Bone punch 2/3/4 any one
9. MI465 Beyor bone rongeur
10. MI477 Lacrimal chisel
11. MI480 Bone guage
12. MI483 Periosteal elevator double ended
13. MI486 Traquair elevator double ended
14. MI489 Mallet s.s.
15. MI492 Nasal speculum
16. MI495 Nasal packing forceps
17. MI501 Pigtail probe double ended
18. MI504 Lacrimal dilator set of 3
19. MI513 Bowman lacrimal probe set of 4
20. MI516 Intubation set with tubing
21. MI612 Jayles fixation cum corneal forceps
22. MI681 Dressing forceps serrated jaws
23. MI687 Moorfield Conjunctival forceps
24. MI825 Westcott tenotomy scissors
25. MI861 Eye scissors delicate
26. MI870 Strabismuss scissors
27. MI918 Needle holder castroviejo long
28. MI933 Kalt needle holder with lock
29. MI951 A.C.Wash cannula
30. MI957 Randolph cyclodialysis cannula
31. MI963 Lacrimal cannula
32. MI1029 Silicon bulb with connector
33. MI1056 Sterilizing case with matts.






Muscle Surgery Set



① ② ③ ④ ⑤

1. MI135 Speculum clerk
2. MI174 Desmarres lid retractor
3. MI210 Castroviejo caliper graduated
4. MI402 Greafe muscle hook
5. MI405 Jameson muscle hook
6. MI423 Cullers muscle hook
7. MI534 Serrafine Bulldog Clamp (2)
8. MI537 Towel forceps
9. MI543 Artery forceps straight
10. MI546 Artery forceps curved
11. MI600 Corneal wide handle toothed forceps
12. MI681 Dressing forceps serrated jaws
13. MI690 Tissue forceps 1x2 teeth
14. MI762 Jamson muscle forceps Right / Left
15. MI765 Worth muscle forceps Right / Left
16. MI825 Westcott tenotomy scissor
17. MI867 Stevens tenotomy scissor
18. MI894 Barraquer needle holder
19. MI951 Wash cannula
20. MI1029 Silicon bulb with connector.
21. MI1050 Sterilizing box with matt.



⑥ ⑦ ⑧ ⑨ ⑩ ⑪





⑫ ⑬ ⑭ ⑮ ⑯



⑰ ⑱ ⑲ ⑳ ㉑





Glaucoma Set




1. MI102 Barraquer wire speculum
2. MI135 Speculum clerk
3. MI210 Castroviejo caliper graduated
4. MI267 Iris hook angled
5. MI345 Patson dissector double ended
6. MI354 Cyclodialysis spatula double ended
7. MI393 Tooks corneal knife
8. MI438 Blade breaker and holder
9. MI471 Kelly Glaucoma punch
10. MI474 Castroviejo holth punch
11. MI498 Trabeculotomy probe right / left
12. MI537 Towel clip cross action
13. MI543 Artery forceps straight
14. MI546 Artery forceps curved
15. MI552 Hess Iris forceps
16. MI588 Girard colibri forceps
17. MI600 Castroviejo corneal forcep
18. MI630 Tying forceps straight
19. MI681 Dressing forceps serrated
20. MI825 Westcott Tenotomy scissor
21. MI843 Vannas scissor angled upward
22. MI855 Barraquer Iris scissors
23. MI924 Castroviejo needle holder
24. MI951 A.C.Wash cannula
25. MI1029 Silicon bulb with connector
26. MI1053 Sterilizing box with matt



Corneal Transplant Set



1. MI132 Guyton Park Speculum
2. MI135 Clerk Speculum
3. MI207 Fleringa ring set
4. MI210 Castroviejo caliper
5. MI222 Castroviejo corneal trephine
6. MI225 Trephine disposable
7. MI237 Corneal radial maker 6 blades
8. MI345 Patson spatula double ended
9. MI390 Paufiguer knife
10. MI534 Serrafine Bulldog clamp (2)
11. MI540 Towel clip
12. MI543 Artery forceps
13. MI552 Hess Iris forceps 1x2 teeth
14. MI561 Hoskin colibri forceps
15. MI573 Suturing forceps
16. MI591 McPherson forceps
17. MI630 Jaffe tying forceps
18. MI792 Westcott tenotomy scissor
19. MI795 Corneal scissors left
20. MI813 Corneal scissors right
21. MI825 Corneal Scissors full curved
22. MI897 Needle holder micro curved
23. MI942 Air Injection cannula
24. MI951 A.C.Wash cannula
25. MI1029 Silicon bulb with connector
26. MI1053 Sterilizing box with matts
27. MI3045 Eye ball stand




Enucleation Set


1. MI135 Clerk eye speculum
2. MI174 Desmarres lid retractor
3. MI336 Bunge Evisceration spoon
4. MI339 Mulle Evisceration spoon
5. MI342 Wells Evisceration spoon
6. MI402 Strabismus hook
7. MI543 Artery forceps (4)
8. MI693 Fixation forceps
9. MI867 Stevens tenotomy scissors straight
10. MI879 Enucleation scissors curved
11. MI1050 Sterilizing box with matt
12. MI3042 Eye collection bottle



Vitreotomy Set




1. MI135 Clerk eye speculum
2. MI162 Localiser
3. MI426 Gas detachment hook
4. MI513 Wilson foreign body forceps
5. MI522 Vitrectomy scissors horizontal
6. MI525 Foreign body forceps vertical
7. MI528 Vitrectomy forceps horizontal
8. MI537 Towel clip cross action (2)
9. MI543 Artery forceps straight
10. MI546 Artery forceps curved
11. MI552 Hess Iris forceps toothed
12. MI564 Hoskins corneal forceps
13. MI600 Corneal wide handle forceps toothed
14. MI615 Kelman McPherson forceps
15. MI708 Superior Rectus forceps
16. MI711 St.Martin corneal cum fixation forceps
17. MI837 Vannas scissor
18. MI903 Needle holder barraquer stoudt
19. MI942 Air Injection cannula
20. MI951 A.C.Wash cannula
21. MI1053 Sterilizing case with matts
22. N.I. Sleeve spreading forceps
23. N.I. Plug holding forceps
24. N.I. Pluggers holding forceps
25. MI375B Membrane scratcher
26. MI372 Membrane peeler
27. N.I.. Pluggers (2)
28. MI372B Membrane spatula, curved
29. MI375A Membrane spatula, pointed tips
30. MI375 Membrane spatula, hooked
31. MI378A Membrane spatula, knob
32. MI1005A Infussion cannula with tubing



Lasik Surgery

1. MI1500 Flap protector & Spatula
2. Mi1503 Flap protector blades angled 45°
3. MI1518 Stabilization ring with swivel handle
4. MI1542 Corneal marker for Lasik surgery
5. MI1548 Lasik spatula & flap unzipper
6. MI1551A Lasik spatula single ended
7. MI1554 Aspirating speculum
8. MI1557 Speculum adjustable spread for lasik
9. MI1560A Buratto lasik flap holding forcep
10. MI1563 Lasik cannula flattened spatulated tip
11. MI1581 Lasik cannula 1 front & 2 side ports
12. MI3081 Marking pen sterilized
13. MD39 PVA Sponge packet of 5 sterilized





Lasek Surgery

1. MI1524 Lasek trephine (8mm with 8.5mm well)
2. MI1524A Lasek trephine (9mm with 9.5mm well)
3. MI1536 Lasek Epi peeler semi-sharp
4. MI1539 Lasek micro hoe & Epi peeler
5. MI1554 Aspirating speculum
6. MI1587 Lasek alcohol cannula
7. MI1530 Epi-Dissector
8. MD39 PVA sponge set of 5 sterilized
9. MI1053A Sterilizing case



friMed

**New
Instruments**



Akahoshi Capsulorhexis Cross Action Forcep
 Cross action mechanism reduces loss of viscoelastic material through the incision/ through side port incision anterior chamber can be maintained and rhexis edge can be better controlled. Capsule can punctured easily and grasp it firmly.

MI-2001



Akahoshi Paddle Prechopper
 Broad blade safely and effectively prechops soft nucleus.

MI-2004



Akahoshi Super Micro Combo Prechopper

- * 20 guage combo prechopper through 1.2mm to 1.5mm incision,
- * Designed to be used through a side port incision, Micro phaco incision and a laser phaco incision,
- * Opening & closing is within the anterior chamber,
- * Reduces phaco time by 50% and for both hard & soft nucleus.

MI-2007



Akahoshi Levitator (for a fallen nucleus)

Helps move nucleus to the bag with holding stable manner, infusion not required.

MI-2010



Akahoshi Ring Levitator

to lift a dropping nucleus from pars plana
 The viscoelastic material can help in levitation.

MI-2013



DSAEK S MARKER

Smark to easily identify the stromal from the endothelial side of the button when inserting and positioning the donor tissue in the recipient eye

MI-2016





New Instruments

DSAEK Endorhexis Forcep
 Designed to remove the Descemet's membrane and endothelium during DSAEK surgery. Blunt tip to prevent stromal damage.
MI-2019



Spatula for Lasik enhancement
 Very easy to use dual instrument for re-cutting, lifting & repositioning the corneal flap during Lasik enhancement.
MI-2022



Micro Scissor
 Squeeze handle
 360° rotatable
MI-2025 Straight
MI-2028 Curved
MI-2031 Vertical



Micro Forceps
 Squeeze handle
 360° rotatable S.Steel
MI-2034 Straight
MI-2037 Curved
MI-2040 Vertical



Sleeve spreading forcep
MI-2043



Pluggers holding forcep
MI-2046



Disposable Range

Disposable Cydodialysis Spatula
MD-03



Disposable Irrigating Vectis
MD-06 Small
MD-09 Large
MD-12 Serrated



Disposable Hydrodissection Cannula
20 to 27 gauge
MD-15



Disposable Hydrolysis Cannula
20 to 27 gauge
MD-18



Disposable Air Injection Cannula
MD-21



Disposable Snare sterilized
MD-24



Disposable Cannula T-Loupe with Silicon-tubing Right / Left
MD-24



Disposable Simcoe Cannula with Silicon tubing
21 / 22 / 23 gauge
MD-27




PVA Instrument Wipe
MD-33



PVA Surgical Points
MD-36




PVA Sponge packet of 5
MD-39





TriMed

Equipments



Fiber optic head light
with cable & light source.
MI-3000



Senior Illuminated Trial Set
Reduced aperture lenses fitted in
metal rims with +/- sph 32 pairs each,
+/- Cyl. 20 pairs each & 10 prisms,
7 accessories consisting 225 lenses
MI-900 - same as above but with
powder coated.
MI-3003



Trial Frame for reduced aperture
Lenses with adjustable P.D. and nose
bridge light weight
MI-3006
MI-3006A - adjustable frame.



Rotter test type drum
with remote control with Duchrome
test, with four dots test for three
meters distance.
MI-3009





Equipments


Internally illuminated rotating drum
distance vision with Duochrome test,
friend list, spotlight, worth four dots test.
MI-3012 - wall mounted manually operated.
MI-3012A - motorised with remote control.



Near vision drum
internally illuminated with stand
MI-3015



Ishihara Charts for colour blindness
MI-3018



Retinoscope Mirror
MI-3021



Cross Cylinder 0.25
MI-3024



Maddox
Wing test for squint
MI-3027



Prism Bar
Horizontal & Vertical Set
MI-3030



Binocular Loupe
Loupe 3.5X with frame
MI-3033



Streak Retinoscope
Operate either with large dry batteries or Rechargeable cell. Separate control for Beam adjustment & Streak rotation 360° Beam Adjustment to Divergent, Parallel & convergent. Large Semi-Reflector interference Mirror. With One Spare Bulb 3 V.
MI-3036



Artificial Eyes,
assorted sizes box of 50
MI-3039



Eye Collection Bottle & Stand
MI-3042
MI-3045





Equipments

Slit Lamp

Technical Specifications : Binocular Microscope
Eye Pieces 10x, 15x , Objective 1x, 1.6x
Total Magnification 10x, 16x and 24x

Illumination : Slit Image Rotation 0° to 180°
Tilting Illumination upto 20°, Filter Disc. Cobalt Filter,
Green Filter, Yellow Filter, Natural density and open aperture
Slit Diaphragm Disc Six apertures of 12, 9, 7, 3 and 0.2mm &
wedge shaped diaphragm for infinitely variable slit lengths.
Halogen Lamp 6 Volt 20 Watts

Standard Accessories : 15x Eyepieces, Replacing Lamp one
piece for illumination unit, Replacing Mirror, Hruby Lens,
Plastic breath shield, Testing Rod, Dust Cover, Fuse
MI-3048



Keratometer

Features : Coincidence focusing system. 15 x eyepieces.
One position instrument. Measure both meridians without
changing optical system. Simple vertical adjustable head to
fit patient. Astigmatic axis 0 to 180 Degree.
The achromatic optical systems to precise measurement of
the corneal radius curvature and accurate measurement of
radius curvature of contact lens. Dual eye-level sighting
system.

Measuring range : Radius of curvature, 9.5 - 6.5 mm
(in 0.25 mm steps) Dioptre 36 - 52 (in 0.25 steps).
Paint finish : Black or White

Standard Accessories : Dust Cover, One spare bulb
Optional Accessories : Manual Stand/motorized Stand
MI-3051



Lensometer

Specifications :
Type Occular type external reading system
Vertex Power Range $\pm 25D$ step : 0.25D step upto $\pm 10D$
0.50D step beyond $\pm 10D$
Cylindrical Axis Scale 0 to 180° (1° step)
Prismatic Power Range 5 in the viewfield (1 setp)
Target Pinshole Corona
Eyepiece Focussing + 3D to -5D
Diameter of Lens Accomodated 30 to 80mm
Tilting Angle 30° to 90°
Light Source AC 220V 3W LED
MI-3054



Doctors Ergonomic chair for O.T. & O.P.D.
Imported D.C. Motor gives as smooth and jerk free
vertical motion during up and down movement.
Touch & stop switch under main top for up movement
of chair gives extra safety.
240V A.C., Total weight of unit (approx.) 250kg.
MI-3057



Electromechanical for Refraction
A strong main arm-pole to accomodate
vision testing equipments, like distance /
near vision drum, projectors, dom light
and additional tray for contact lens.
Provision for slit lamp, keratometer,
trial lens set etc.
MI-3060





Equipments

Wetfield Surgical Bipolar coagulator cautery supplied with one insulated forcep, crazer silicon cable & handle with foot switch.

MI-3063



Super Pinky for reducing preoperative intraocular pressure.

MI-3066



Operating Microscope

Microscope Section

Optical Head Indigenous,

Eye Piece Wide field 10x individually adjustable

Eye Inclination 45° (ophthalmic) 90° (ENT & general)

Interpupillary Distance Adjustable from 54mm to 75mm

Working Distance 200mm, **Total Magnification** 5x to 25x variable

Field of Vision 12mm to 52.5mm

Illumination Section Type

Fiber cable cold light guide 24V 250 Watts incorporated with heat filter.

Field of Illumination 55mm

Stand Section Type Floor stand, Movable with castor wheels,

Vertical movement of counter balance arm 300mm

Fine Focussing Manual 40mm movement or motorised with foot pedal

Radius of rotation arm 360° **Base Size** 570mm - 590mm

Weight 45Kg. approx

OPTIONAL ACCESSORIES Beam Splitter with CCD Camera attachment, CCD Camera, Table Clamp attachment.

MI-3069



Synophthophore

Movement of Optical Tubes

Horizontal : Adduction + 50° -- Adduction -40°

Vertical : Hyper 30° -- Hyper 30°

Torsional : Incyclo 20° -- Excyclo 20°

Slide Illumination Rheostat controlled 12V Lamp for each slide.

After Image illumination by 12V Lamp (for better illumination)


Auto Flashing Auto Flashing of Slide illumination either Simultaneous or Alternate in Rapid & Variable modes.

Mode & Mode Selection Five Modes of Slide illumination namely Normal, Flashing Right, Flashing Left, Flashing R+I & Autoflashing, can be selected by a single selector knob.

Standard Accessories A set of slides containing 9 Pairs power Cord. Dust Cover. Two Spare Bulb 12V.

Optional Accessories Instrument Table (Motorized)/ Manual Table Haidinger Brushes

MI-3072



Binocular Indirect Ophthalmoscope
MI-3075



SCHIOTZ Tonometer
MI-3078



Marking Pen Sterilized
MI-3081



Instrument Tray
MI-3084



Kidney Tray
MI-3087



Bowl
MI-3090






Equipments

Sterilizer
MI-3093



IRIS Retractor sterilize
MI-3096



Tension Injector
MI-3099



Capsular Tension Ring
MI-3102



Back Flush. PSD
MI-3105



Silicon Tires
MI-3108



Silicon Sponge Rods
MI-3111



Silicon Bands
MI-3114

MI-3114

Corneal Loupe
10x folding type
MI-3117



Punch Trepine
MI-3120




Spectacle Magnifier
MI-3123



Universal Eye Shield
Green Colour autoclavable
MI-3126



Eye Curved Needle
Pkt of 6
MI-3129



Cheatele forcep
MI-3132



Sponge Holding Forcep 8"
MI-3135



Autoclave Vertical
MI-3138





TriMed Technologies

Customise Surgical Solution & Instruments

4/E-1, Second Floor, Nuchana Tower,
Alaknanda Shopping Complex, Alaknanda, New Delhi - 110019 (INDIA)

Ph:- +91-11-41051188, Mobile:-+91-8459130431

Email: trimedtechnologies@gmail.com, info@trimed.co.in

Web: www.trimed.co.in