

DRYSTAR 5302

Can reliability and affordability really fit together?



High-quality images

everywhere

and anywhere

When you want a day-to-day hardcopy imager that meets all your imaging needs, look no further than the DRYSTAR 5302. Its robust design supports long life and reliable operation, while the two on-line trays and table-top size provide optimal flexibility and usability. It fits where you want it and adapts to your imaging needs, providing the affordable and reliable answer to your CR/DR printing requirements. Plus, you get the all high-value DRYSTAR family features, including Agfa's A#Sharp-enhanced, award-winning Direct Digital Imaging technology (DDI); ecological and easy operation and maintenance; and straightforward, DICOM-native integration. DRYSTAR 5302 is the flexible solution to your daily tasks.

Fits where you need it


DRYSTAR 5302 offers high-quality printing in a convenient table-top size, so you can offer next-to application printing in even the smallest spaces. This freedom of placement, combined with the low investment and running costs, makes DRYSTAR 5302 the perfect match for many modalities, and the right answer for your printing needs. You make the choice: use it dedicated to one modality, or easily connect it to a small, multiple-modality network. While it is primarily intended to support CR/DR, it can also be used with CT/MR or you can use it as a replacement printer for an existing wet imaging solution.

Flexibility for maximum usability


DRYSTAR 5302's two adaptable media trays can handle all film types and sizes, and both are available to you, on-line, all the time, for great operational efficiency and flexibility. The trays handle blue base and clear base (Agfa dry film) media in 5 different sizes, for remarkable adaptability (8 x 10 inch, 10 x 12 inch, 11 x 14 inch, 14 x 14 inch and 14 x 17 inch). All in all, the DRYSTAR 5302 is the printer that conforms to your needs, not the other way around.

Consistency and quality: a dynamic pair

DRYSTAR 5302, like all the Agfa family of DRYSTAR imagers, incorporates Direct Digital Imaging technology. With its solid-state technology, it provides controlled, consistent image quality in a stable, long-lasting platform. Using DRYSTAR 5302 with Agfa's state-of-the-art dedicated media guarantees diagnostic quality for all applications.



Agfa's DDI has been enhanced to include A#Sharp technology, for even sharper image quality. Equipped with this powerful combination, the DRYSTAR 5302 provides you with diagnostic comfort and quality, ease of use, reliability and long life, consistently. DDI and A#Sharp: a matching set.



Networked or dedicated... it always fits in



Friendly to the environment and your workflow

DRYSTAR 5302's DDI technology is fast, clean and environmentally-friendly. DDI's simple mechanism means less wear and tear on your equipment. There is no chemistry, no processors, no darkroom and no expensive chemical disposal costs. Total daylight media handling during loading or maintenance increases the ease of use and comfort while reducing costs.

The affordability answer

Using DDI, DRYSTAR 5302 has no complex or expensive moving optical parts. You can take full advantage of the benefits of trouble-free upkeep and maintenance. From the initial investment, through the low maintenance and repair costs, the DRYSTAR 5302 offers a remarkably low total cost of ownership. Less hassle, less trouble, less cost: that's a simple choice!

Seamless integration into the heart of your facility

DRYSTAR 5302 has been developed to fit right into your facility's system. The DICOM-native imager makes network connectivity easy.



For more information on Agfa, please visit our website on www.agfa.com ■

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert NV, Belgium, or its affiliates. DRYSTAR is a trademark of Agfa-Gevaert NV, Belgium, or its affiliates. All rights reserved. All information contained herein is intended for guidance purposes only, and characteristics of the products and services described in this publication can be changed at any time without notice. Products and services may not be available for your local area. Please contact your local sales representative for availability information. Agfa-Gevaert NV diligently strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

© 2020 Agfa NV
All rights reserved
Published by Agfa NV
Septestraat 27 - 2640 Mortsel
Belgium

NGLNH GB 00202001

