



TriMed Technologies

Customise Surgical Solution & Instruments

4/E-1, Second Floor, Nuchana Tower,

Alaknanda Shopping Complex, Alaknanda, New Delhi - 110019 (INDIA)

Ph:- +91-11-41051188, Mobile:-+91-8459130431

Email: trimedtechnologies@gmail.com, info@trimed.co.in


Web: www.trimed.co.in

Expand your vision.

The **Murdoch Mechanical Arm** has the stability and strength that are required to hold the Nathanson Liver Retractor in situ. The arm can be resterilized by using a steam autoclave and is supplied nonsterile.

Order Number	Reference Part Number	Description
G03364	MMA-1000	Murdoch Mechanical Arm
Replacement Parts		
G49849	MMA-1000-BHS	bolt and handle set
G38698	MMA-1000-MBC	main body clamp (engraved C)
G38700	MMA-1000-MCW	main collar washer (engraved B)

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Murdoch
MECHANICAL ARM

Nathanson Liver Retractors are specially designed to fully support and retract the liver during laparoscopic upper GI surgery in order to improve the view of the operating field.

The shape of the stainless steel device evenly distributes the weight load, and the anodized steel eliminates reflections. The retractors come in two sizes and can be resterilized by steam autoclave.

Order Number	Reference Part Number	Description	Distal tip distance mm
G03123	NLRS-1001	large	75
G03395	NLRS-1002	extra large	95

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Nathanson

LIVER RETRACTOR



TriMed Technologies

Customise Surgical Solution & Instruments

4/E-1, Second Floor, Nuchana Tower,

Alaknanda Shopping Complex, Alaknanda, New Delhi - 110019 (INDIA)

Ph:- +91-11-41051188, Mobile:-+91-8459130431

Email: trimedtechnologies@gmail.com, info@trimed.co.in

Web: www.trimed.co.in